



THE GREAT OUTDOORS

Permobil's X850 has rugged good looks and is designed for rough terrain. Featuring rear-wheel drive and fully independent suspension, the X850 allows you to climb obstacles up to four and a half inches high, ideal for the outdoors. Its directional headlights turn with you, and the spring-loaded rearview mirrors ensure you do not get caught on obstacles along the trail. Offering servo-steering with a joystick, the X850 may be programmed to meet the user's requirements with a 9-mph maximum speed. Plus, the Corpus[®] seating system is the world's most popular powered seating system, and customization is almost limitless with depth, width, height, and angle of multiple seat components being easily adjustable.

X850
CORPUS

POWER FUNCTIONS



Power adjustable tilt

SHROUD COLORS



Vintage Wine



Crystal Blue



Black Onyx



Key protected control box



Independent rear suspension



Optional 180° rotational seat for easy transfer from either side. Available in power or manual options.



Power extendable front end allows custom fit for every user and added stability.



Power adjustable recline

SPECIFICATIONS

Maximum Speed	9 MPH	Gradient Capacity (uphill and downhill)	15°
Driving Range	22 - 28 Miles	Side Slope	12°
Max User Weight	300 lbs	Power Seat Elevator	8"
Ground Clearance	3.25"	Seat Widths	17" - 21" (Adjustable)
Battery Capacity	Group 31 97 Amps/Hr	Seat Depths	14" - 22" (Adjustable)
Width	28"	Backrest Heights	21" - 29.5"
Length	49" - 57"	Seat To Floor Height	22"
		Electronics (standard)	SafeGate



THE POWER OF MOBILITY

Permobil Inc.
 6961 Eastgate Boulevard
 Lebanon, Tennessee 37090
 800.736.0925
 www.permobil.com